



## DeepUnity Viewer for Pathology

UNCOVER THE FULL  
POTENTIAL OF A DIGITAL  
WORKFLOW IN PATHOLOGY

DeepUnity Viewer for Pathology is a module  
of the medical device DeepUnity Viewer



Improve collaboration  
within your institution

### Viewing Capability



Multidisciplinary  
by design



Digitalized workflow  
for more efficiency


DIAGNOSTIC IMAGING IT

## Pathology viewing capabilities – now part of our Universal Viewer

DeepUnity Viewer for Pathology supports  
pathologists in their daily work routine.  
Pathologists can have quick access to  
digitalized slides, analyze them and share all  
digitalized medical content in an easy way  
within their institution.



## Archiving Capability



Fast, robust and scalable archiving



Open standard communication (FHIR/HL7) and multi-tiered storage



Automatically DICOMization of slides



## DIAGNOSTIC IMAGING IT Pathology archiving capabilities

For a fully comprehensive solution, we offer our archive solution DeepUnity Medical Archive on top. Providing a robust and scalable solution for your needs. Right after the slides are scanned, they are converted into DICOM, enriched with patient information from the Laboratory Information System, and stored in DICOM.